MANUEL D' INSTRUCTIONS

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Models	CQ-12	CQ-27	CQ-29	CQ-33	CQ-36
Style	desk model	floor model with five-point base and casters	Same as left model	Same as left model	Same as left model
Diameter of mineral plate	124 mm	124 mm	166 mm	124 mm	166 mm
Extendible range	16" (40 cm) vertical, 10" (25 cm) horizontal adjustment with 2 sections, spring loaded arms	24" (60 cm) vertical, 24" (60 cm) horizontal adjustment with 2 sections, spring loaded arms	Same as left model	Same as left model	Same as left model
Emission headlamp direction	Any direction within a half sphere may be tilted up to 90 degrees. Device headlamp may be rotated 50 degrees	Any direction within a half sphere may be tilted up to 90 degrees. Device headlamp may be rotated 360 degrees	Same as left model	Same as left model	Same as left model
Telescoping main post	fixed	allows 12" (30 cm) vertical adjustment	Same as left model	Same as left model	Same as left model
Distances to floor	6" (15 cm) minimum to 18" (45 cm) maximum	12" (30 cm) minimum to 60" (152 cm) maximum	Same as left model	Same as left model	Same as left model
Treatment distance	8" to 12" (20 cm to 30 cm)	8" to 12" (20 cm to 30 cm)	Same as left model	Same as left model	Same as left model
Weight	12 lbs (5.5 kgs.)	25 lbs (11 kgs.)	Same as left model	Same as left model	Same as left model

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Symbols	Meaning
	Manufacturer
	Authorized Representative in the European Community
	Symbol for the marking of electrical and electronics devices according to Directive 2002/96/EC. The device, accessories and the packaging have to be disposed of waste correctly at the end of the usage. Please follow Local Ordinances or Regulations for disposal.
	CE 0123 mark of conformity The device is complied with the Essential Requirements for Medical Devices according to EC Directive 93/42/EEC
	Serial number
	Refer to instruction manual
	Warning: Hot surface
	Protective earth (ground)
	Year of Manufacture
	ON (Power)
	OFF (Power)

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Guidance and manufacturer's declaration - electromagnetic immunity			
The CQ-12 is intended for use in the electromagnetic environment specified below. The customer or the user of CQ-12 should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge(ESD) IEC 61000-4-2	± 6 KV contact ± 8 KV air	± 6 KV contact ± 8 KV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 KV for power supply lines	± 2 KV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 KV differential mode ± 2 KV common mode	± 1 KV differential mode ± 2 KV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the CQ-12 requires continued operation during power mains interruptions, it is recommended that the CQ-12 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Mains power quality should be that of a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

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Guidance and manufacturer's declaration – electromagnetic immunity			
The CQ-12 is intended for use in the electromagnetic environment specified below. The customer or the user of CQ-12 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3V _{rms} 150kHz to 80 MHz 3V/m 80 MHz to 2.5 GHz	3V _{rms} 3V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the CQ-12 including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CQ-12 is used exceeds the applicable RF compliance level above, the CQ-12 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the CQ-12.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

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For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer.

NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Guidance and manufacturer's declaration -electromagnetic emission		
The CQ-27/29/33/36 is intended for use in the electromagnetic environment specified below. The customer of the user of the CQ-27/29/33/36 should assure that it is used in such and environment.		
Emission test	Compliance	Electromagnetic environment -guidance
RF emissions CISPR 11	Group 1	The CQ-27/29/33/36 use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The CQ-27/29/33/36 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

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Guidance and manufacturer's declaration-electromagnetic immunity			
The CQ-27/29/33/36 is intended for use in the electromagnetic environment specified below. The customer or the user of CQ-27/29/33/36 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge(ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood,concrete or ceramic tile.If floor are covered with synthetic material,the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2kV common mode	±1 kV differential mode ±2kV common mode	Mains power quality should be that f a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60%dip in U_T) for 5 cycles 70% U_T (30%dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60%dip in U_T) for 5 cycles 70% U_T (30%dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the CQ-27/29/33/36 requires continued operation during power mains interruptions, it is recommended that the CQ-27/29/33/36 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Mains power quality should be that of a typical commercial or hospital environment
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

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Guidance and manufacturer's declaration – electromagnetic immunity			
The CQ-27/29/33/36 is intended for use in the electromagnetic environment specified below. The customer or the user of CQ-27/29/33/36 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the CQ-27/29/33/36, including cables, than the recommended separation distance calculated from the equation applicable
Conducted RF IEC 61000-4-6	3V _{rms} 150 kHz to 80 MHz	3V _{rms}	<p>to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer and d is the recommended separation distance in metres(m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3V/m 80 MHz to 2.5 GHz	3V/m	
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>^a Field strengths from fixed transmitters, such as base stations for radio(cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CQ-27/29/33/36 is used exceeds the applicable RF compliance level above, the CQ-27/29/33/36 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the CQ-27/29/33/36.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

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Recommended separation distances between portable and mobile RF communications equipment and the CQ-27/29/33/36			
<p>The CQ-27-29/33/36 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled, The customer or the user of the CQ-27/29/33/36 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CQ-27-/29/33/36 as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter(m)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	$d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			